

## COX4I1-cytochrome c oxidase subunit IV isoform 1 Gene

Reference: PTXBC008704

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX IV-1; COX4; COX4-1; COXIV; COXIV-1; cytochrome c oxidase subunit 4 isoform 1, mitochondrial; cytochrome c oxidase polypeptide IV; cytochrome c oxidase subunit IV; cytochrome c oxidase subunit 4I1		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	cytochrome c oxidase subunit IV isoform 1
<b>GENE ID</b>	1327
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC008704
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	COX4I1
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon )

atgttgctaccaggtatttagcctagttgcaagcgagcaattccacctctgtgtgtacgagctcatgaaagtgttgaagagcgaagactttcctcccagcttatatgatcggcgtgaccaccttgcggaggtggcccatgcaagcacctgtctgccagcagaaggcaatgaaggagaaggagaagcctctctggagcagcctcctcatggatgagaaagtcagttgtatcgattaaagttcaaggagactttgctgagatgaacagggcctcgaacgagtggaagacggtgtggcgtgccatgttctcatcggttcaccgcgctcttatcatgtggcagaa  
gcactatgtgtacggcccctccgcaaaccttgacaaagagtggtggccaagcagaccaagagatgctggacatgaaggtgaacccatcaggccttagcctccaagtgaggactacgaaaagaacgagtggaagaagtga