



Reference: PTXBC008725

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VATI; synaptic vesicle membrane protein VAT-1 homolog; membrane protein of cholinergic synaptic vesicles; vesicle amine transport protein 1 homolog (T. californica); vesicle amine transport 1		

For research use only.

GENE DESCRIPTION	vesicle amine transport protein 1 homolog (T. californica)
GENE ID	10493
ORIGIN SPECIES	Human
ACCESSIONS	BC008725
VECTOR	pDONR223
NCBI SYMBOL	VAT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]