

ANXA8L2-annexin A8-like 2 Gene

Reference: PTXBC008813

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANXA8L2; ANXA8; VAC-beta; bA145E20.2; annexin A8-like protein 1; annexin A8-like protein 2; annexin A8L2; annexin VIII; annexin-8; vascular anticoagulant-beta; annexin A8-like 1		

For research use only.

GENE DESCRIPTION	annexin A8-like 2
GENE ID	244
ORIGIN SPECIES	Human
ACCESSIONS	BC008813
VECTOR	pDONR223
NCBI SYMBOL	ANXA8L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgggtgaaagcctggattgaacagagggtgtcacagtgaagagcagctccactcaaccagaccctgatgcagagaccctctacaagccatgaaggggatcgggtgcgggtccaactgctcagccaccaagcagctgccttcgctccctcctccgcctcaccag
tgtgtcaccctggggcagcaggggtcactgtgtgttaacctgcagggaacagcagcaggtatcatgatgtgtcctaccaagagaagcaacacgcagcggcagcagatcgcaagtcctcaaggctcagttcggcaaggacctcactgagacctgaagtctgagctcagtgcca
agtttgagaggctcattgtggcccttatgtatccgcatagatacagaagccaaggagctgcatgacgccaatgaagggcagcagggatgatgtgagcagcttgggacccggcactggccctccaagacgcacaggaatctgtatcgccgagggcgaagatctgtgggactgatgag
atgaaattcatcaccatcctgtgcacgcgcagtgccactcacctgtgagagtgaaatgcacccaaaacctccacagctactttgcagagagactctactatgcatgaaggagcagggacgcgtgatgggacctgataagaacatcgtttcaaggagcagattgacttaaatctatc
aaatgtcacttcaagaagatgtacggcaagacctcagcagcatgatcatggaagaccagcggcgactacaagaacgccctgctgagcctgggtggcagcgaccctga