

AKR1B10-aldo-keto reductase family 1, member B10 (aldose reductase) Gene

Reference: PTXBC008837

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AKR1B11; AKR1B12; ALDRLn; ARL-1; ARL1; HIS; HSI; aldo-keto reductase family 1 member B10; ARP; SI reductase; aldo-keto reductase family 1, member B10 (aldose reductase); aldo-keto reductase family 1, member B11 (aldose reductase-like); aldose reductase-li		

For research use only.

GENE DESCRIPTION	aldo-keto reductase family 1, member B10 (aldose reductase)
GENE ID	57016
ORIGIN SPECIES	Human
ACCESSIONS	BC008837
VECTOR	pDONR223
NCBI SYMBOL	AKR1B10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacgtttgtgagctcagtaccaagccaagatgccattgtggcctgggcaacttggagctctctcttggcaagtgaagaagcagtgagggtggccattgatgcagatcgccacattgactgtgcctatgtctatcagaatgaacatgaagtgggggaagccatcaag
agaagatccaagagaagcgtgtgaagcgggaggacctgttcacgtcagcaagttgtggccaactttcttgagagacccctgtgaggaaagcctttgagaagaccctcaaggaccctgaagctgagctatctggacgtctatcttacttggccacagggttcaagtctgggatgac
ctttcccaagatgataaagtaatacctcgttggaaaagcaacgttcttgatgcctggaggccatggaggagctggtgatgaggggctggtgaaagcccttgggtctccaattcagccacttccagatcgagaagctcttgaacaaactggactgaatataaaccagtga
taaccaggttgatgtcaccatacctcacacaggagaactgatccagtactccactccaaggcacccttacggcctacagccccctgggctctccgtagacattggccaagccagaagacccttccctgctggaggatccaagattaaggagattgtcgaagcaca
aaaaccgcagcccagggttctgatccgtttccatccagagggaatgtgattgtcatcccaagctgtgacaccagcagcagctgttgagaacattcaggtctttgactttaattgagtgatgaggagatggcaaccatactcagcttcaacagaactggaggcctgtaacgtgttgcaat
cctctcatttgggaagactatccctcgtatgcagaatattga