

SLC25A44-solute carrier family 25, member 44 Gene

Reference: PTXBC008843

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	solute carrier family 25 member 44		

For research use only.

GENE DESCRIPTION	solute carrier family 25, member 44
GENE ID	9673
ORIGIN SPECIES	Human
ACCESSIONS	BC008843
VECTOR	pDONR223
NCBI SYMBOL	SLC25A44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggagacaacgcaacatcagatcatcagatgggaacacctggacaagaagaagtctactgtttggtgtggcaatgacaatgatgccgtgtcagtgctaccattaccctcatccgacccggttcaagttcagaagggaagagcctctaccatggaccttcgatgcct
catcaagatcctcgagcagatggtatcactggccttaccgaggggttcctggtcaataccctaccctcatctctggccagtgttatgtcaccactatgagctcaccggaagttttagctgactacagccagagtaaacagctcaaatcactggtggctggtgctcagcctccctgtgg
cccagagcatcacagtgccattgatgtatctcccagaccctgatgatgcaacgcaagggtgagaaaaatggccgtttcagggtgcggggaacccagagggacaagggttagtgccttggccaaccaaaggacatcatcagcagatcctgagcctgatggactcgcggcttc
taccgagcctatgtgcttcaactatcccaaacagtgctgtctgtggcccttctactctatgagagcagctctctacctgtctcctaggagtgccctcactgtcttcaagctgtctcgggcccctggctgagcactgcctccatctcacaatccatggatgtc
atacgaaccctgtgaggttgaggcaagaactccatcctcagcctcagacagctgatggcagaagaaggccttggccctcatgaaggccctctcggccagaatctcagccacacctccaccattgtcattgtgtggctgatgagagcctcaagaactcagcctccgac
ctgagctgtggactcagacactggtaa
```