



Reference: PTXBC008847

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-dependent RNA helicase DDX24; DEAD (Asp-Glu-Ala-Asp) box helicase 24; DEAD (Asp-Glu-Ala-Asp) box polypeptide 24; DEAD box protein 24; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 24; S. cerevisiae CHL1-like helicase; DEAD-box helicase 24		

For research use only.

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 24
GENE ID	57062
ORIGIN SPECIES	Human
ACCESSIONS	BC008847
VECTOR	pDONR223
NCBI SYMBOL	DDX24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]