

NOP10-NOP10 ribonucleoprotein homolog (yeast) Gene

Reference: PTXBC008886

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOP10 ribonucleoprotein; snoRNP protein NOP10; NOP10 ribonucleoprotein homolog; DKCB1; NOLA3; NOP10P; H/ACA ribonucleoprotein complex subunit 3; homolog of yeast Nop10p; nucleolar protein 10; nucleolar protein family A member 3; nucleolar protein family A		

For research use only.

GENE DESCRIPTION	NOP10 ribonucleoprotein homolog (yeast)
GENE ID	55505
ORIGIN SPECIES	Human
ACCESSIONS	BC008886
VECTOR	pDONR223
NCBI SYMBOL	NOP10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttccagctattacctaagcagcaggagatcgagctctatacgctgaagaaattgacccgatgggacaacagacgtgcagccatcctgctcggttctccagatgacaataactctcgacaccgaatcaccatcaagaacgcttcaagggtgctcatgacccagcaaccgcccgtgctctga