

SEDLP-spondyloepiphyseal dysplasia, late, pseudogene Gene

Reference: PTXBC008889

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	spondyloepiphyseal dysplasia, late, pseudogene
GENE ID	10597
ORIGIN SPECIES	Human
ACCESSIONS	BC008889
VECTOR	pDONR223
NCBI SYMBOL	SEDLP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggaagcttctactttgttaattgttggtcaccatgataaccagttttgaaatggagttttgccagctgggaaggcagaatccaagacccatcgtcatctgaaccagttcatgctcatgctctcgacctcgtagatgagaacatgtggctgtcgaacaacatgtacttgaaaact
gtggacaagtcaacagtggtttgtgtcagcattgtcaccgcgggcatatgaggtttattatgcttcacataagacaagaagatggaataaagaactctttactgatgtttatgattatataaaaatttcaatgaatccattttatgaaccaattctcctattcgatcaagtcatttgaca
gaaaagttcagtttcttgggaagaacaccttttaagctga