

## ZBTB38-zinc finger and BTB domain containing 38 Gene

Reference: PTXBC008901

|                |                                                                                                                                                                   |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CIBZ; PPP1R171; ZNF921; zinc finger and BTB domain-containing protein 38; protein phosphatase 1, regulatory subunit 171; zinc finger and BTB domain containing 38 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | zinc finger and BTB domain containing 38                                                          |
| <b>GENE ID</b>          | 253461                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC008901                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | ZBTB38                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgcaggaggagccttgcacaggggaatgacccagaaccagtgagagacgcccactcggccttgcacatccgagtgcatggagatgagtgagtgctgacgcaagtaccaggattccactgacaaacgtggcgcccttactacaactacaaacccaaaaaagaatcca  
gacagttgaaaaaaatgaggaaagtcaactggaggagagcagcggaaacaggagcccgagccataaatgtaataccagcagaactggattgcgcctggggaaggctcctcaggataaaccttggaggaagaactaaagagatcccaagctgcagtgtaactctgtg  
atggagacaaagcagtgggggctggaaccaagggaagggccacacatcttactctcggccatagctcgcagctctgcgccaagcagttccagagccctccacactcaaatgcacatgagatgtcacaccggggagaagccataccagtgcaagacctgcggacggtgctt  
ttcggfcaaggaaacttacagaacatgaacgcatccacctgggcttgaaggagttctgtcagtgattgcaacaaggcattcaccttgatgagaccctcaaatccatgaagaatccatactggagaaaagcgttaccactgtcagttctgtttcagagattttgtatctccacaaa  
aggaaatcacgagcagagcatattcgggagcataatgggaaggctatgcttccagtgcccaaaattgcacaaacagctgctgcttggaatgcacaaaagaacacttatcaaaagcccaagtcagcaggagaaaatagggtgacgtgtgccacgaaactcaaatcccttg  
gagaatcaacatttattggttcagaagacaatgaccaaaggataacatacaaacgggtgtggaaaatgtgtccttga