

ARF6-ADP-ribosylation factor 6 Gene

Reference: PTXBC008918

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADP-ribosylation factor 6; ADP ribosylation factor 6		

For research use only.

GENE DESCRIPTION	ADP-ribosylation factor 6
GENE ID	382
ORIGIN SPECIES	Human
ACCESSIONS	BC008918
VECTOR	pDONR223
NCBI SYMBOL	ARF6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggggaaggctctcaaaatctcgggaacaaggaaatcgggatcctcatgtggcctggacgcggccgcaagacaacaatcctgtacaagtgaagctggccagtcggtgaccaccattcccactgtgggttcaactggagacggtgacttcaaaaatgtcaagttcaacg
tatgggatgtggcggccagacaagatccggcctctcggcgttactactacactgggaccaaggctcctcctcgtagtgactgcgccgaccgcgaccgcatcgcagctcgcaggagctgcaccgattatcaatgaccgggatgaggacgccataatcctcctc
gccaacaagcaggacctgccgatccatgaaaccccacgagatccaggagaactggcctgacccgattcgggacaggaactgtatgtcagccctcctgtccacctcaggggacggactctatgagggctcatggttaaccttaactacaaatcttaa
```