

LOC389833-similar to hypothetical protein MGC27019 Gene

Reference: PTXBC008921

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to hypothetical protein MGC27019
GENE ID	389833
ORIGIN SPECIES	Human
ACCESSIONS	BC008921
VECTOR	pDONR223
NCBI SYMBOL	LOC389833
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctgccctagagccgtttctgtaggacattcatggcgtttccaggacctctgccagctccattggagttttgggtccctccacccagtaggtgcatacttctgccctccagccagctcacaccccccaacaccccaagctggtgagtgaggagctg
tacaagtcctcacagcgtgctgaggagaagcttacaagcggagcagctccctgcgcacctcaaggacatccacatgagcctggagaaggacgtcactgccatgaccaacagctcttctcatcgaccgccagaagtcagtgcccatgtactgtctacccaccattctgcagctggt
ggctaccagtga