



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RABAC1-Rab acceptor 1 (prenylated) Gene

Reference: PTXBC008950

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

PRAF1; YIP3; prenylated Rab acceptor protein 1; PRA1 domain family 1; PRA1

family protein 1; Rab acceptor 1 (prenylated); prenylated Rab acceptor 1; Rab acceptor 1

For research use only.

GENE DESCRIPTION	Rab acceptor 1 (prenylated)
GENE ID	10567
ORIGIN SPECIES	Human
ACCESSIONS	BC008950
VECTOR	pDONR223
NCBI SYMBOL	RABAC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcgcagaaaggaccagcagaagaatgccgagcgcggaaggctgacgcgccacgacctctgtccgaagctgattccctccggcgagcgccggagtgctgagagcggcgccgcgaccatccggccctggagcaccctctgggaccagcagcgtcttcaaggccccc
caactctgggagagctgtgccaagcgcctcgtacgcaacgtggagctactacagagaacaattgtgtctgtctgggacctatcctgactgtgtggcgacgtcccatattgtgtgctggctctggctgtcttttggcgctgttaacctctatctctgcgcacctggagatccaagcttgctg
ctttggccgagaggtgagcccgccatcatgatgtcctgctggcgaggcactctctcccctctctctgctggctgctgctggcgctcggcgctgctgggtgctggagaccacctgggtgctatcgctgccacgctgcttccaccagattaggctgtgacggggagagagctgca
gcttccttcctgctgta

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1