

C16orf45-chromosome 16 open reading frame 45 Gene

Reference: PTXBC008967

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C16orf45; MINP; migration inhibitory protein; chromosome 16 open reading frame 45		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 45
GENE ID	89927
ORIGIN SPECIES	Human
ACCESSIONS	BC008967
VECTOR	pDONR223
NCBI SYMBOL	C16orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaattaaagcaatcttctccaccatctggaagccgagagcctctgaggcgctatggggcgggtggaggagacggcttgaaaacggagagactggggagaaatcagctggacatcatccatggcgagagacaacatgatccagaggagattgagctggagatggcaaaa
attcagcgtctcgggaagtcttggtccggggagctctgagctcaggttcattgatggatgacatccagctctgcaaggacatcatggacttgaagcaggagctgcagaacttggtccatccagaaaaagaaaaaaccaactgcagaagcagagagagatgagctaatccagaa
gatccacaactggtgcagaagagagacttctggtgacgatcgagggtcgagcgggttaaggagagcaagaagaagacaaggaaatggctgatttcctgagaatcaagttaaacctctagacaaagtaaccaaatctccagccagctccgggcagagaagaagcagagccccc
acctagcaagccacggtggccaagacggggctggcactgatcaaggattgttcggggccaccagtcgaacatcatgtag