

YRDC-yrdC domain containing (E. coli) Gene

Reference: PTXBC008984

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	yrdC N6-threonylcarbamoyltransferase domain containing; yrdC domain containing; yrdC N(6)-threonylcarbamoyltransferase domain containing; yrdC domain-containing protein, mitochondrial; DRIP3; IRIP; SUA5; dopamine receptor interacting protein 3; hIRIP; isc		

For research use only.

GENE DESCRIPTION	yrdC domain containing (E. coli)
GENE ID	79693
ORIGIN SPECIES	Human
ACCESSIONS	BC008984
VECTOR	pDONR223
NCBI SYMBOL	YRDC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacgctcggaggagctcaacaaggacctaaccctttacgcctctgttagcattcgattcctgatcatgctttatgcaagactggctcagatgtttgagggctcgttctctcactagtccaacctcagctccaggccagtctctgaatgtcaggagttccaggatctctggtcctagttgtccttggttattgatggggacaaattgggatggccagagccccgagtgctgccttggtcctcaactgtggttgattgtctgtgccggaaagtgttgcatcattgtccaggctgtgccctggaaagtactacagccatcctccaacagaagtacggactgctccctcacatgcgtcctacgtgta