

C9orf30-chromosome 9 open reading frame 30 Gene

Reference: PTXBC008993

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome 9 open reading frame 30
GENE ID	91283
ORIGIN SPECIES	Human
ACCESSIONS	BC008993
VECTOR	pDONR223
NCBI SYMBOL	C9orf30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaaacacgaattataaagcctgccaatactctcagaattggaaaagacatctgctgctttagtagaaaagtataaatatgtgctggaatgtaagaaaagtgatgcgcgaactattgcccttaagcagcgtacctggcaggcgtgcccacgaatacaactctcagccagc
gtgtccctgcgggatttcaaacagctgaagaagtgtctgggagaacatcaaggctcggacaaaaaattatggcccatgaaaaggagagaaaagtgaacggagcgtcagccctctcctgagtagccacgtcctagggaaggagaagatgccagcatgctgccggagcagctctact
tctgcagagccccccggagggagcccgataccaccccgacgctcagcccaagaatcattgtctgtttcaatagagaactgtgcgatgatgagaagaggttcatacattttcagtagtgaggggacctctcaacctgaacctcgtgttcagctgtcagaataacagccaataaa
aactacaggagcaaacctctcagggaaggtgttttaaaaaagatgcatgagggaagaacacccaataaagtgccatcttcaactgcaactgatacaaatgaatgagggtgcatgtggccaaaatccagcagatagagcgagagtgtagatggcagagggaacacaggataaaaat
ggaagtctcaataaaaagaagatgtattgggaaagaaaactacaacttttaccaaggatggcctgtttcctcatttaaccggcccttccaattcgccctaa