

## FAM164A-family with sequence similarity 164, member A Gene

Reference: PTXBC009074

|                |                                                                                                                                                                              |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | protein FAM164A; FAM164A; C8orf70; CGI-62; zinc finger C2HC domain-containing protein 1A; family with sequence similarity 164, member A; zinc finger C2HC-type containing 1A |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | family with sequence similarity 164, member A                                                     |
| <b>GENE ID</b>          | 51101                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC009074                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FAM164A                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggaggagctggagagagatggaggtgtgtccagttggagaattgtaccttgcagatttgggaagaacattctccagtagcattaaaaaacatggaccatttgcagaagactgcaactaaaaacggaagactttgattcaagcagacagagagctgaaggaaactgatattc  
caacagtaaaacgctcaaacaggccagaccaccaagaacacatcaattggagaaggaaacatgaagaattcattgctaccataagagcagctaaaggccttgatcaggccctcaaaagggtggcaaaacttctctctctcctcactcttctatgatctgattatattcaatgtcc  
atattgtcagaggagattcaatgaaatgcagctgatagacataataattctgtaagaacaggcagcagctatttagtaataaaggaaatttctgcagataccaaaaggaaacaaacttctcgacacaggtgtataagccaccgcaacttaaaagtcaaatctctggaactgcatcat  
caggatctcagattaccgcagccaaagtggcgtggcaaaactgtttaggtgtcttcaggtaagtgtcttcaagtagcagctcttgggaacaacactcagacattgtctccctctcataaaggatagcagccctcatgcaggagctaatgtcaaacccgaaattccacaccact  
agtttggcaagaatcctgcccaggtgtgttacaacaaaagaaaacataatactgagagctacatagccagcagatggggactgtgcatctcccttaattggtggaataattaagcattgaaggacattcacttgaaacttaccaaaattctgcatgagtggtggactaaatacc  
ctgtagaattggccaaatttctgtgaatgtggcattcgaagaatgattctatga