

NDUFA13-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 13 Gene

Reference: PTXBC009189

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B16.6; CDA016; CGI-39; GRIM-19; GRIM19; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 13; CI-B16.6; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 13; NADH-ubiquinone oxidoreductase B16.6 subunit; cell death regulatory protein GRIM-19; c		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 13
GENE ID	51079
ORIGIN SPECIES	Human
ACCESSIONS	BC009189
VECTOR	pDONR223
NCBI SYMBOL	NDUFA13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgtaaaaggtgaagcaggacatgcctccgcgggggctatggcccatcgactacaaacggaactgcgcgtcgaggactgtcgggctacagcatgctggccatagggttggaacctgatctacgggactggagcataatgaagtggaaacctgagcgcaggcg
cctacaaatcgaggacttcgaggtcgatcgcgctgttgcactgttacaggcagaaacgacgggagacgttcagatgcttcgggagaaacctggaggaggagccatcatgaaaggacgtgccgactggaaggtggggagctgtgtttccacacacccgtgggtgcc
cccttgatcggggagctgtacgggctgcaccacagaggaggtctcctcatgccagccacggcttcattgtgtacacgtag