

MRPL13-mitochondrial ribosomal protein L13 Gene

Reference: PTXBC009190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L13; L13A; L13mt; RPL13; RPML13; 39S ribosomal protein L13, mitochondrial; MRP-L13; mitochondrial ribosomal protein L13		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L13
GENE ID	28998
ORIGIN SPECIES	Human
ACCESSIONS	BC009190
VECTOR	pDONR223
NCBI SYMBOL	MRPL13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgagtttcttagggcgccccagcaatggccacttttctagaatatgtatctcttagatgggaaatgcagccacctggcaacttgcctatggcactataagactcagggattacataaacctgtgtaccatgcactgagtactgtgggatcatgtgttataatgaacacaa
gacacattgcattttctggaacaaatgggaacaaaagtatactcttcgcatactggctaccaggtggatttagacaagtaacagctgctcagcttcacctgagggatccagtggaattgtaaaactagctattatggcatgctgccccaaaaccttcacagaagaacaatgatggaaag
gttgcattctttccagatgagtatattccagaagatattcttaagaatttagtagaggagcttctcaaccacgaaaaatacctaacgtctagatgagtacacacagaagaataagacgcttcccaagattgtggactccacctgaagattatcggtataa