



Reference: PTXBC009194

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HP55; MTPOLB; PEOA4; POLB; POLG-BETA; POLGB; DNA polymerase subunit gamma-2, mitochondrial; DNA polymerase gamma accessory 55 kDa subunit; mitochondrial DNA polymerase subunit gamma-2; mitochondrial DNA polymerase, accessory subunit; p55; polymerase (DNA		

GENE DESCRIPTION	polymerase (DNA directed), gamma 2, accessory subunit
GENE ID	11232
ORIGIN SPECIES	Human
ACCESSIONS	BC009194
VECTOR	pDONR223
NCBI SYMBOL	POLG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1