

## C5orf30-chromosome 5 open reading frame 30 Gene

Reference: PTXBC009203

|                |                                                                                             |                     |                |
|----------------|---------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | UPF0684 protein C5orf30; UNC119-binding protein C5orf30; chromosome 5 open reading frame 30 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 5 open reading frame 30                                                                |
| <b>GENE ID</b>          | 90355                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC009203                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C5orf30                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggaagtcgatattatggagagctagaagtagccctgaccaccttgccttcctgggctgaggccaactccccgggaaaggcggaggcagagaagccccgctgctccagcacacctgtccccgatgcggaggaccgtgtcagggtaccagatcctacacatggacttaacta  
tttgggttggcttcacgactggcgaggaactcctgaagtagctcagaagtgacagagggtgaagagagcaaaagcagaagccatgccatccttacgctccaaacagctagatgcaggacttgcctctctgtttgtataaaaccagaagtaggtactaccagccatcagagattccag  
ctgtcaatggcaggaggcgaaggcggatgccaagctcaggagacaagtgactaaatctttaccttgaaccttacaaggccctccatgggctctgctcttctttaaaggtaagagggtcactccaatctctggactacctcaatctagataaatgatcaaggagccagctga  
tacagaagtgtacagtaccagcttcaacacctaaccctccgaggggaccgtgtgtttgtaggaataatacatga