

CTD-2267G17.3-X antigen family, member 2-like Gene

Reference: PTXBC009232

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	X antigen family, member 2-like
GENE ID	728242
ORIGIN SPECIES	Human
ACCESSIONS	BC009232
VECTOR	pDONR223
NCBI SYMBOL	CTD-2267G17.3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggcgaggagatcaacatagcctagcccaagaagaagtttacagcctcctgagctgattgggctatgctgaaccactgatgaagagcctaagaagagaaccaccctaaaagtcggaatcctacacctgacagaagagagaagatgatcagggtgcagctga
gattcaagtgcctgacctggaagccgatctccaggagctatgcagacaaagactgggagtgatggaagggtggaagctgaaggggaagattctacaaaagcagagcattttaaatgccagaagcagggtgaagggaaatcacagggttaa