

## FDFT1-farnesyl-diphosphate farnesyltransferase 1 Gene

Reference: PTXBC009251

|                |                                                                                                                                                    |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | DGPT; ERG9; SQS; FPP:FPP farnesyltransferase; presqualene-di-diphosphate synthase; squalene synthetase; farnesyl-diphosphate farnesyltransferase 1 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | farnesyl-diphosphate farnesyltransferase 1                                                       |
| <b>GENE ID</b>          | 2222                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC009251                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | FDFT1                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggagttcgtgaaatgcctcgccacccgaagagttctacaacctggtgcgttcggatcggggcaagcgggaaggtgatgccagatggaccagactcgtcagcagcagcctgaaaactgtctacaggtatctcaatcagaccagtcgagtttcgagctgttatccaggcg  
ctggatggggaatgcgcaacgcagtgatgcatatttatctggtctccgagctctggacacactggaagatgacatgacctagtggtgaaagaaggtcccgctgttacacaactttcactcttcttaccacagactggcggttcagtgagagcaaggagaaggatgccagggtg  
ctggaggacttccaacgactcctctgagtttagaaatctggctgagaataccaacagtgattgccgacatttccggagaatgggcatgggatggcagagtttttgataagcatgtgacctctgaacaggagtgaggacaagtactgccactatgttgcctggctggcgggaattggcc  
ttcccgcttttctcagcctcagagtttgaagacccttagttggtgaagatacagaacgtgccaaactctatgggctgtttctgcagaaacaaacatcatccgtgactatctggaagaccagcaaggagggaagagagttctggcctcaagagtttggagcaggtatgttaagaagttaggg  
gattttcgaagccggagaatattgacttggcctgcatgctgcaatgaactataaccaatgcactgcaccacatccagatgcatcacctaccttgcagactcagaaccagagtggttttaactctgtgctatccacaggtgatggccattgccacttggctgcctgttataataacca  
gcaggtgttcaaaaggcagtgaaattcggaaaggcaagcagtgacctgatgatggatgccaccaatatgccagctgtcaaaagccatcatatcagtatatggaagagatttatcatagaatccccgactcagaccatctctagcaaaacaaaggcagatcatccaccatccgga  
cgcaaatctccaactgtcagctgattccgaagccactactccccatctactgtctgttgcacgttttggctgcctgagctggcagtaactgaccactctctccaggtaacagaagactatgttcagactggagaacactga