

C17orf59-chromosome 17 open reading frame 59 Gene

Reference: PTXBC009261

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf59; PRO2472; BLOC-1-related complex subunit 6; lysosome-dispersing protein; lyspersin; BLOC-1 related complex subunit 6		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 59
GENE ID	54785
ORIGIN SPECIES	Human
ACCESSIONS	BC009261
VECTOR	pDONR223
NCBI SYMBOL	C17orf59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgtgagcagcaggagggaagacaacgacgaggaggcgccggcgagcagagccggcgctgttctccagccgtctcaggacagccgagcctggacgggctgagcaggcgtgctggcgccgggtctcaggagtgccgagtcgagcggcgccggcgccgacccatctccagtcctcggagctcgaaggcacagtgagccgacggcgaccttaccactttgtcgcacaacctgcaactcaagatccgtctgagcggcgccctccaccccgcttctgcccctgagcggccctgccagcgcctgcacccacaccaca
ccggccattcccccatcgaccccgaggtgctgaggatctgagcgggtgagtcgggagctgggagccgggtggacgctctctcgggtctggtggcgccggtgcaggagctgacggcgctgagcgtgggtgacccagacctaccgcgatgctgtggactccttaggtgaa
gccgtggacatgagcatcaaggcatgtacacccctgctggcgctgctgaggagctggagcgggctctgcagccgggtcaggggctggcgccgaagtcgggatalccgactactctgaggtgttggagccctgtgcaagtga