



Reference: PTXBC009282

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AHC2; ATP1A1; CAPOS; DYT12; RDP; sodium/potassium-transporting ATPase subunit alpha-3; ATPase, Na+/K+ transporting, alpha 3 polypeptide; Na(+)/K(+) ATPase alpha(III) subunit; Na(+)/K(+) ATPase alpha-3 subunit; Na+, K+ activated adenosine triphosphatase al		

[illegible]