



## Antibodies and Proteins

## ProteoGenix SAS

## Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

### CD7-CD7 molecule Gene

Reference: PTXBC009293

|                |                                                                                                                                                  |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                              | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                             | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CD7 molecule; CD7 antigen (p41); T-cell antigen CD7; GP40; LEU-9; TP41; Tp40; T-cell leukemia antigen; T-cell surface antigen Leu-9; p41 protein |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | CD7 molecule                                                                                      |
| <b>GENE ID</b>          | 924                                                                                               |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC009293                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CD7                                                                                               |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

atggcgcggcctcgcaggctcgtcgtcgtccctcgtcttcgtcgtggctcgcgcgcctgcctggggccctggctgcaccaagagtgtagcagctctcccactgcacgactgccccgtgggagctccgtcaacatcacctgctccaccagcggggcctcgtcggggtacctga  
ggcagctcggggccacagcccaagacatcattactacgagacgggggtgggtcccactacgacagacgggttcggggccgctcgtactctcagggctccaggacaacctgactatcacatgcaccgctgcagctgtcggagactggacactcacctgccggccatcacgg  
aggctcaatgtctacggctccggacactcgtcgtcgtgacagaggaacagctccaaaggatggacaaatgtcgcgcgcgcccaagggcctcgtccctccctgcccccagcagaggctcgcgcctcctgaccgcgacagacgtcctgcctcctgaccgccagacagctcctgcctcctgaccgccagcagcctcgt  
ccctccctgcgcgcctcgcgggtgatctcctctcctcctcgggctgggctgggggtggcgtgtgtcgtcgtgcgacacacagataaagaacctgtcgtcgtgcgcgggataaagtgcgcgcgcatgtgtgtgtacgaggacatgtgcacagccgctgcaaacagctgtctctcccccac  
caggcaggatga