

UBB-ubiquitin B Gene

Reference: PTXBC009301

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-50; polyubiquitin-B; epididymis secretory protein Li 50; polyubiquitin B; ubiquitin B		

For research use only.

GENE DESCRIPTION	ubiquitin B
GENE ID	7314
ORIGIN SPECIES	Human
ACCESSIONS	BC009301
VECTOR	pDONR223
NCBI SYMBOL	UBB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagatctctgtgaaacccctaccggcaagaccatcaccttgagggtggagccagtgacaccatcgaaatgtgaaggccaagatccaggataagggaagcattccccgaccagcagaggtcatctttgaggcaagcagctggaagatggccgtactcttctgactacaac
atccagaaggagtcgaccctgcacctggtctgctgctgagaggtggtatgcagatctctgtgaagacctgaccggcaagaccatcaccttggaagtggagccagtgacaccatcgaaatgtgaaggccaagatccaggataaagaaggcatccctcccgaccagcagaggtca
tctttgcaggcaagcagctggaagatggccgacctcttctgactacaacatccagaaggagtcgacctgacctggtctgctgctgagaggtggtatgcagatctctgtgaagacctgaccggcaagaccatcactctggaggtggagccagtgacaccatcgaaatgtgaagg
ccaagatccaagataaagaaggcatccccccgaccagcagaggtcatctttgaggcaagcagctggaagatggccgaccttctgactacaacatccagaagagtcgacctgcacctggtctgctgacctgaggggtggctgttaa