

TMEM129-transmembrane protein 129 Gene

Reference: PTXBC009331

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D4S2561E; E3 ubiquitin-protein ligase TM129; transmembrane protein 129, E3 ubiquitin protein ligase; transmembrane protein 129		

For research use only.

GENE DESCRIPTION	transmembrane protein 129
GENE ID	92305
ORIGIN SPECIES	Human
ACCESSIONS	BC009331
VECTOR	pDONR223
NCBI SYMBOL	TMEM129
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagcccgaggtagccttcaactcgcctatctgggtgttcgcggtgtgttcacgccaacgagttccacgcgccggggctcacggtagcagaacctgtgtcgggctggctgggcagcaggagcgcgccttcgtgccctccacttgcgcgcacggcccaacgtg
ttgtgccactcgtgtgcccgtcggtactactgtggcatgtgcttgcggcttcagaaaagcgctccacgcccagcagggccctgaggcctggcgctcttctgtgtggtgacccctccctccatgcctgcatcctgactactactgtccctgacgggtggcctg
ccacccactggcgccacccctggccctctacgccctccacagctgtggctggcaggtgtgttcctctctgtcaacactgagttccggcggttgacaagtttgccaccggtagcaccaggtgcccgtgtgattgtgacagacacgtgggtgatgaaggtaaccactaccgagtgacgtg
gcccagcagcaggagcgtgacctgactgtgacggagtctcggcagcatgagctctgcccagactcgaactggccgtgcagctcctcaccatccgtgtggccagcaccacccctgtgtgagccttgacatctgagctggaggcctgcatag