

C21orf70-chromosome 21 open reading frame 70 Gene

Reference: PTXBC009341

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf70; PRED56; protein FAM207A; family with sequence similarity 207 member A		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 70
GENE ID	85395
ORIGIN SPECIES	Human
ACCESSIONS	BC009341
VECTOR	pDONR223
NCBI SYMBOL	C21orf70
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaagtgggggttcgcgcgccgagtcaccaggtgcccgtgaggcgaaaggaggcccccggcccgccccctgccccggaggcgaccctcgcggcctcgccgcccggggaaggactggcgcttcacacaccaacatctttgccaggaccaagat
agacccagcgccttgggtcagaagctggagctggagctgaggtgctactccgtcaggagaggtgaggcaggtcagtgacggagcgtcccttcacaggagaggtgagaggccaagaccgttttggccaaaggagaaaatgaagctgaggcgtgagcaatgggttc
agaaaatcgaagccataaaactggctgagcagaagcacaggaggagcgaggcgaggccacgggtggtggtggggacctgcaccctcagggatgccctgcccagctgctgggctcagggctggcagccggcgccaaagccgcagcaggagagcaacaagcccc
ggccctcagagctcagcggatgagcgcagcccagagacagcagcttcgagggaagaaggaccgggttcaggagctgctggccagtcggcctacagagcagccccctgggtggccatcgggcagacgtggccccgcagatgcagctggaagatggcgccagctctga