

## CDCA7L-cell division cycle associated 7-like Gene

Reference: PTXBC009352

|                |                                                                                                                             |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | JPO2; RAM2; cell division cycle-associated 7-like protein; transcription factor RAM2; cell division cycle associated 7 like |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cell division cycle associated 7-like                                                            |
| <b>GENE ID</b>          | 55536                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC009352                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | CDCA7L                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggagttggcgaactcgtaccagatccctaaagaagtggctgacatcttaacgccccagtgatgatgaagattgttggctccgagatgatgtccatggaaccctctctgcagaggagagctgcgatatgttactactagatcagggaacacagcaggatgctgctttcattc  
caaatcttcacagaagagctaaagaattttatagaggacactgactcagagactgaggattttgcaggattttacgagagtgatctgaatggaaagactaacccagaagtaattggtcgtggagtcagatttgatgatgagcgaagcatctttggtgagcagaggagaggaagatg  
aagaagaagataaggctacccctagaagaagcaggtctagaagaagtagtattggtcttcgagtagcctttcagttccccaccaagaagctggccacaacaacagataaaaaacagttctccgagcagttgtttctagcgcacgcttacagaatgagaaaaaacaattcttgaagaaga  
aagactgtagacaggtgatacaaaagggaagattctacctctgagctcaggatgactctcgggatgagagccaggagagttcagatgcttctgctgaaaaaggaccatgaacatcaaggagaacaagccatgcttgcctcatttggcggaattgaactcgatgccagatttctccagta  
cgaaccccaacctcagctctaggaagaagacagtgaggcggttccttcggaggacagatcacgcggcgtatgaacccaacccggagtgccggtcctctgagaagtttgcctagagaacttcactgtctcagcgttaatttgcgggaagattttacagcttcgaagaaggag  
acaattgggggaaatgccggagtagacagacgacgtcaccttatcttcttcggccagtgaggatataccgaagaggacttagaaaatgttgcataactgttcgagataaaatctatgataaagtttgcgttaacacgtgccatcaggttcgacaaaagaccatcgacaccaagaca  
gtgtgtcgggaaccagggttctgtgtgtgtgcgaggacagttctgtggacatgcctgcggaaaccgctatgggagggatgtcagatcgccattgctggaccggattgggtgtgtccccctgtcgtgggatctgaattgcagctactgtcggagcgtgacgcccgtgtgccaagga  
atctctcattctgtgccaagtttatggttatgacaattgtaaggaatatctggagagcttacaagaaggagctggtagaagacaataa