

TAGLN2-transgeline 2 Gene

Reference: PTXBC009357

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HA1756; transgeline-2; SM22-alpha homolog; epididymis tissue protein Li 7e; transgeline 2		

For research use only.

GENE DESCRIPTION	transgeline 2
GENE ID	8407
ORIGIN SPECIES	Human
ACCESSIONS	BC009357
VECTOR	pDONR223
NCBI SYMBOL	TAGLN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaacaggggacatgcatatggcctgagccgggaggtgcagcagaagattgagaacaatatgatcagatctggagcagatcctgatccagtggatcaccacccagtgccgaaggatgtggccggccccagcctggacgcgagaactccagaactggctcaaggatgg
cacggtgctatgtgagctcattaatgacactgtaccccgaggggcagggccccagtaagaagatccagggcctccacatggccttcaagcagatggagcagatctctcagttcctgcaagcagctgagcgtatggcattaaccactgacatcttccaaactgtggacctgtggaagga
aagaacatggcctgtgtgcagcggacgtgatgaatctgggtggcgtgagccgagatgatgggctcttctctggggatcccaactggttccttaagaatccaaggagaatcctcggaacttctcgataaccagctgcaagagggcaagaacgtgatcgggttacagatggg
caccaaccggggcgctctcagggcagcatgactgctacgggatgccacgccagatcctga