

TOMM22-translocase of outer mitochondrial membrane 22 homolog (yeast) Gene

Reference: PTXBC009363

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1C9-2; MST065; MSTP065; TOM22; mitochondrial import receptor subunit TOM22 homolog; mitochondrial import receptor Tom22; translocase of outer membrane 22 kDa subunit homolog; translocase of outer mitochondrial membrane 22 homolog; translocase of outer mit		

For research use only.

GENE DESCRIPTION	translocase of outer mitochondrial membrane 22 homolog (yeast)
GENE ID	56993
ORIGIN SPECIES	Human
ACCESSIONS	BC009363
VECTOR	pDONR223
NCBI SYMBOL	TOMM22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgccgcgtgctgctgccgtgcagggaacccagtcctccggacgaattgctccgaaaggcgacgcggagagcctgaggaggagctggaggagcagcagatgaggagctagatgagaccctgtcggagagactatggggcctgacggagatgttccggagagg
gtccggtccgcggcggagccacttttgatcttccctctttgtggtcagaaaatgtacagggtttccaggcgagccttggattgggaccacttcttattgatcctgttcttccgtgttcttggagacggagaagttgcaaatggagcaacgacgcaactgcagcagcggcagatactt
ctaggacctaacacagggctctcaggaggaatgccaggggctctaccctcacttctctggaagatctag