



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DDX47-DEAD (Asp-Glu-Ala-Asp) box polypeptide 47 Gene

Reference: PTXBC009379

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E4-DBP; HQ0256; MSTP162; RRP3; DEAD (Asp-Glu-Ala-Asp) box polypeptide 47; DEAD box polypeptide 47; DEAD box protein 47; E4-DEAD box protein; DEAD-box helicase 47		

For research use only.

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 47
GENE ID	51202
ORIGIN SPECIES	Human
ACCESSIONS	BC009379
VECTOR	pDONR223
NCBI SYMBOL	DDX47
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcaattcttggcccttcgaaaaaaaccacataataataatgcaactctctggctgactgattgaccacttggaaaatacgaaggtttcaactctgagagctctcaaatacttggctgatggaagccgaccgaatactgaatatgatttggagacagaggttgacaagatctcctcaagtg
attcttcgagatcggaataacattctctctctctgcacacatgacgaagcaaggtttcaaaaactcagcgagacagctctgaagaatctctgtaaatgttgcctgttctcttaataatcacagacatggttgaataatcagaacatataattttatctctcttaaatcaaggataactctctgttataattc
aaatgaatctgcctggaaactcttataatattctgcgaatactgataataataccgagacagctgttctctatcgctgaaacttgcctctcatcgccatccctccatggcacaataatgactgagataagcgcctgagataatgaaatggaacccgcctctatcttctatgcaact
gacgttggccagcggaggttggacatacctctatgatagtgtgtgttgcacttgcagacttctaccatctccacatccaaggataatcatcatcctgagtatgtctgaacagctagacgtggcgctccggaagcgctattactttttgcacacagtatgatgtggaactcttccagcgcatagaacacttaattg
ggagaataactacaggttcttccacacaggaatgagatgatttgcgtctgacgacgaacgcgtcgtctgaagcccaaggtttgccgaatggaggttaagggaagcatggagaaaaagaagaacgctcgcgagagagctgtggagataatgactacacagagaggtgtctattgtgtcaggaa
caagtgctgcctgcgaataaatagaacacgcgaagagccctgttaa

More info: Please contact us for any question: contact@proteogenix.fr