

RPS17-ribosomal protein S17 Gene

Reference: PTXBC009407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA4; RPS17L; RPS17L1; RPS17L2; S17; 40S ribosomal protein S17; ribosomal protein S17		

For research use only.

GENE DESCRIPTION	ribosomal protein S17
GENE ID	6218
ORIGIN SPECIES	Human
ACCESSIONS	BC009407
VECTOR	pDONR223
NCBI SYMBOL	RPS17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggccgcttcgacacaaaccgtgaagaaggcggccgggtcatcatagaaaagtactacacgcctgggcaacgactccacacgaacaagcgcgtgtgcgaggagatcgccattatccccagcaaaaagctccgcaacaagatagcaggttatgtcacgcatctgatgaag
cgaattcagagaggccagtaagaggtatctccatcaagctgcaggaggaggagagagaaaggagagacaattatgttcctgaggtctcagccttgatcaggagattatgaagtagatcctgacactaaggaaatctgaagcttttgacttcggcagctgtgtcaaccttcaggctact
cagcctacagttgggatgaatttcaaaacgcctcggggacctgtttga