



Reference: PTXBC009409

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	truncated TACSTD2; EGP-1; EGP1; GA733-1; GA7331; GP50; M1S1; TROP2; tumor-associated calcium signal transducer 2; 40kD glycoprotein, identified by monoclonal antibody GA733; cell surface glycoprotein TROP2; cell surface glycoprotein Trop-2; epithelial gly		

For research use only.

GENE DESCRIPTION	tumor-associated calcium signal transducer 2
GENE ID	4070
ORIGIN SPECIES	Human
ACCESSIONS	BC009409
VECTOR	pDONR223
NCBI SYMBOL	TACSTD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]