

ZDHC8-zinc finger, DHHC-type containing 8 Gene

Reference: PTXBC009442

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DHHC8; ZDHC1; ZNF378; membrane-associated DHHC8 zinc finger protein; zinc finger protein 378; zinc finger, DHHC domain like containing 1; zinc finger DHHC-type containing 8		

For research use only.

GENE DESCRIPTION	zinc finger, DHHC-type containing 8
GENE ID	29801
ORIGIN SPECIES	Human
ACCESSIONS	BC009442
VECTOR	pDONR223
NCBI SYMBOL	ZDHC8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagaggcaccagggccccaccgtctctgacacagcctgtgggtccccgaccgagtggtccccgccaggctactcctaactaacggttgctttcacggaccccgctggaagcttgtag