

AIF1-allograft inflammatory factor 1 Gene

Reference: PTXBC009474

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AIF-1; IBA1; IRT-1; IRT1; allograft inflammatory factor 1; interferon gamma responsive transcript; ionized calcium-binding adapter molecule 1; protein G1		

For research use only.

GENE DESCRIPTION	allograft inflammatory factor 1
GENE ID	199
ORIGIN SPECIES	Human
ACCESSIONS	BC009474
VECTOR	pDONR223
NCBI SYMBOL	AIF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccaaaccagggtttacagggaggaaaagcttcggactgctgaagggccagcaggagagaggctggatgagatcaacaagcaattcctagacgccaaatagcagtgatgaggatctgccctcaaaactggaaggcttcaagagaaatcatggagttgaccttaat
ggaaatggcgatattgatatcatgtccctgaaacgaatgctggagaaactggagtgcccaagactcacctagagctaaagaaatgaattggagaggtgtccagtggtccggggagagcttcagctaccctgactttctcaggatgatgctgggcaagagatctgccatcctaaaaatgatc
ctgatgatgagggaaaagcgagagaaaaggaaaagccaagggcccccagcagaagaagctatctctgagttgccctga