

MRPL33-mitochondrial ribosomal protein L33 Gene

Reference: PTXBC009475

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf1; L33mt; MRP-L33; RPL33L; 39S ribosomal protein L33, mitochondrial; mitochondrial ribosomal protein L33		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L33
GENE ID	9553
ORIGIN SPECIES	Human
ACCESSIONS	BC009475
VECTOR	pDONR223
NCBI SYMBOL	MRPL33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcctcgcgcgtctcttcccaagagcaagtcacaaacattctggtgagaatggtgagcgaagctgggacaggtttctgctcaacaccaagagaaaccgactgcgggaaaaactgactctttgcattatgatccagttgtgaacaaagagtcctctctgtgaaaagaaaaaat
acgctcccttaa