

DCUN1D1-DCN1, defective in cullin neddylation 1, domain containing 1 (S. cerevisiae)
Gene

Reference: PTXBC009478

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DCNL1; DCUN1L1; RP42; SCCRO; SCRO; Tes3; DCN1-like protein 1; DCN1, defective in cullin neddylation 1, domain containing 1; DCUN1 domain-containing protein 1; RP42 homolog; defective in cullin neddylation protein 1-like protein 1; squamous cell carcinoma-		

For research use only.

GENE DESCRIPTION	DCN1, defective in cullin neddylation 1, domain containing 1 (S. cerevisiae)
GENE ID	54165
ORIGIN SPECIES	Human
ACCESSIONS	BC009478
VECTOR	pDONR223
NCBI SYMBOL	DCUN1D1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacaagttgaatcatgcagaaggataaagtcgtcagtttatgatcttcacacaatctagtgaacacagcagtaagttgtcttctcaaatgactggaagttagatgttgcacagataatctttccaaatcctgaactttatatacagagagagtgtaaaaggatcattggacaggaa
gaagttagaacagctgtacaatagatacaagaccctcaagatgagaataaattggaatagatggcatacagcagttctgtgatgacctggcactgcagccagcattagttgtgattgcatggaagttcagagcagcaacacagtcgcagttctcacaacaggagttcatggat
ggcatgacagaattaggtgacagcatagaaaaactaaaggccagatacccaagatggaacaagaattgaaagaaccaggacgatttaaggattttaccagtttacttttaatttgcaagaatccaggacaaaaaggattagatctagaatggccattgcctactggaacttagtgc
taatggaagatttaattcttagacttatggaataaattttgttgaacatcataaacgatcaataacaaaagacacttgaatcttcttttagacttcagtagatgattgcagatgacatgtctaattatgatgaagaaggagcatggcctgttctattgatgactttgggaattgcacgcctcaa
attgctgggacaaaaagtacaacagtag