

TTL3-tubulin tyrosine ligase-like family, member 3 Gene

Reference: PTXBC009479

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tubulin monoglycylase TTL3; HOTTL; tubulin tyrosine ligase-like family, member 3; tubulin--tyrosine ligase-like protein 3; tubulin tyrosine ligase like 3		

For research use only.

GENE DESCRIPTION	tubulin tyrosine ligase-like family, member 3
GENE ID	26140
ORIGIN SPECIES	Human
ACCESSIONS	BC009479
VECTOR	pDONR223
NCBI SYMBOL	TTL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagggtatcgaacatctggatcgtgaagccaggagccaagtcccgaggagcatcatgtgcatggaccactggaggagatctgaagctgggaacggcaacccgtggatgaaggacggcaagtgggtggtgcagaagtattgagcggccctcctcatctt
ggcacaagtgtgacctcagacagtgtgtctgtgtaactgactggaaccacttacgtgtgtgtctacgcgacagctatatccgtttccacgcagccctctccctgaagaacctggacaagtga