

SEC61G-Sec61 gamma subunit Gene

Reference: PTXBC009480

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SSS1; protein transport protein Sec61 subunit gamma; Sec61 gamma subunit; protein transport protein SEC61 gamma subunit; Sec61 translocon gamma subunit		

For research use only.

GENE DESCRIPTION	Sec61 gamma subunit
GENE ID	23480
ORIGIN SPECIES	Human
ACCESSIONS	BC009480
VECTOR	pDONR223
NCBI SYMBOL	SEC61G
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcaggtaatgcagttgttgagccaagtcggcagttgttaaaggactccatcggtgttaaaagatgcactaaacctgatagaaaagaattccagaagattgccatggcaacagcaataggatttgctataatgggattcattggctctttgtgaaattgatccatattcctattaataac
atcattgttggtgctga