

C21orf56-chromosome 21 open reading frame 56 Gene

Reference: PTXBC009497

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf56; speriolin-like protein; spermatogenesis and centriole-associated protein 1-like protein; spermatogenesis and centriole associated 1-like		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 56
GENE ID	84221
ORIGIN SPECIES	Human
ACCESSIONS	BC009497
VECTOR	pDONR223
NCBI SYMBOL	C21orf56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtccggcctaagaaggtgtgtttctcggagagcagcctcccacccgggacaggaccaggaggagctactacatcaatgagatccagagcttcgcccggcgcgagaaggacgcgcgtgtgtggcgagatgccttcagctggaccgccatcctggcctacgtgttccc
gggcgtgacgcggctctacggttcacgggtggccaacatccccgagaagatcgagcagacctccaccaagtctctggacggctccgtggagcagaggaagctgcgcgagctgacgcagcgtactctggcctgagcgcgcgtggagaagctgggtacagccgcgacgtgca
ccggcggttcagcagttctcatcaacacctacggaatcctgaagcagcggcccgacgtgcgcgcaacccctgcacagcagccggccgcgtgcgcaagctggtcatcgacgtgtgtgcccccaagttcctggcgactcgtgctgctcaactgcctgtgcgagctctcc
aaggaggacggcgaagccctcttcgctgtga