

MRPL20-mitochondrial ribosomal protein L20 Gene

Reference: PTXBC009515

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L20mt; MRP-L20; 39S ribosomal protein L20, mitochondrial; mitochondrial ribosomal protein L20		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L20
GENE ID	55052
ORIGIN SPECIES	Human
ACCESSIONS	BC009515
VECTOR	pDONR223
NCBI SYMBOL	MRPL20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcttctcaccgcgagctctggctgcggaatcgcgtcaccgaccgctactttcggatccaggaggtgctgaagcacgcagcagcacttcggggaaggaaaaatcgtctacaggttgccggtcagaaccgtgattcgagcctttgtgaaatgcacaaagcccgatacctgaag
aaaaagaacatgaggaccctctggattaatcgaattacagctgctagccaggaaacatggactgaagtatccagcgtcattgggaatttagttaagtgccaggtggagctcaacaggaaagtcctagcggatctggccatctacgagcaaaagactttcaaatctttgctgccttgccagt
aggaggcgacacgaaggatttgcctgctgccttgggggatgggaaggaaacctgaaggcattttccagagtggtgcagtaccactga