

CALML4-calmodulin-like 4 Gene

Reference: PTXBC009516

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NY-BR-20; calmodulin-like protein 4; serologically defined breast cancer antigen NY-BR-20; calmodulin like 4		

For research use only.

GENE DESCRIPTION	calmodulin-like 4
GENE ID	91860
ORIGIN SPECIES	Human
ACCESSIONS	BC009516
VECTOR	pDONR223
NCBI SYMBOL	CALML4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggccatgaggtgcctggggccagcccagcggagggtgcagcggcacctgcagaccacgggtagacggaaatggagagctggatttcctcactttctgaccattatgcacatgcaataaaacaagaagacccaagaagaattcttagccatgttgatgttggacaaggagaagaaggttacgtcatggcgtccgacgtgcgtcaaaactcacgagctcgggggagaagctcaccacaaggagtgatgatctctcagggaagcagatgcgaaccaatggcaagtgaagtatgatgaatttatccacaagatcaccttctggacgggactattga