

CXorf40B-chromosome X open reading frame 40B Gene

Reference: PTXBC009523

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein CXorf40B; chromosome X open reading frame 40B		
For research use only.			

GENE DESCRIPTION	chromosome X open reading frame 40B
GENE ID	541578
ORIGIN SPECIES	Human
ACCESSIONS	BC009523
VECTOR	pDONR223
NCBI SYMBOL	CXorf40B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtttgctgcctcctccggcagccttatgctggcctttgtttaaataagaaatgaatcaagactgtggagacgcgtggcgtcctctgctgagcagccagcgggaactgtaccatcgccgtccacattgctcacagggaactgggaaggcgatgcctgtcgggagctgctggtggagagact
cgggatgactcctgctcagattcaggccttgcaggaagggaagtttggcagggatgatacgggactcgtgacattggggaacttgcattgccccgaagacttaactcccgatgaggttgggaactagaaatcaagctgcactgaccaacctgaagcagaagtacctg
actgtgatfttcaaacccagggtggttactggagccatacctaggaaaggagcgaaggatgtattccaggtagacatccagagcactgatcccttggggcatgaagtgtga