

CDADC1-cytidine and dCMP deaminase domain containing 1 Gene

Reference: PTXBC009562

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NYD-SP15; bA103J18.1; cytidine and dCMP deaminase domain-containing protein 1; protein kinase NYD-SP15; testis development protein NYD-SP15; cytidine and dCMP deaminase domain containing 1		

For research use only.

GENE DESCRIPTION	cytidine and dCMP deaminase domain containing 1
GENE ID	81602
ORIGIN SPECIES	Human
ACCESSIONS	BC009562
VECTOR	pDONR223
NCBI SYMBOL	CDADC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagaagctgggcagatgcaaaatctggagagcgcgagggccggcggtcagtcagcaccagactggcagcatgaccggtcagataccaaggcttttaagtcaaccttttcaactctgctcagcctctggatggagctcttccagcagaagccagcggcaaaatctcagaaatgaagagggaagcatggaccttaggagataatgaagagagagaccagagatctctctgacaaaagacagaaaacatgttctgctgtttgaaatgattgtaaatgctggagttaa