

SNRNP35-small nuclear ribonucleoprotein 35kDa (U11/U12) Gene

Reference: PTXBC009622

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HM-1; U1SNRNPBP; U11/U12 small nuclear ribonucleoprotein 35 kDa protein; U1 snRNP-binding protein homolog; U11/U12 snRNP 35 kDa protein; U11/U12 snRNP 35K; U11/U12-35K; protein HM-1; small nuclear ribonucleoprotein 35kDa (U11/U12); small nuclear ribonucle		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein 35kDa (U11/U12)
GENE ID	11066
ORIGIN SPECIES	Human
ACCESSIONS	BC009622
VECTOR	pDONR223
NCBI SYMBOL	SNRNP35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaccagctaacatgaacgattggatgccatgccaggagatgatccactcaaaagcgggcagcattgatggcaccgatgaagaccacacgaccgcgcggtctggaggccaatgctggcacgatatgtcccaacaaagggtgcataggagatccctcctcacctgtttgtg
gccagactaaacttcagaccaaggaggacaaatfaaaggagctctttcccgctatggtagacatccggcgcttcggctggtcagggaacttggtcacagggttttcaaaagggtacgccttcacgaatacaaggaggagcgtgccgtgatcaaaagcttaccgagatgctgatggcctggt
tatggaccagcatgagatatttggactacgagctggaaggactctcaaaagggtggatccctcggcgacttggagcggctcttgggggaaaaaaggagctctggcaactgagatttggggacgggacccgctttcgaacattataacttggcagttgttaaaacgacctata
gagagggaacgggaaaggcgggagcgcgtctcgtcccgagaagacactgggactcgaggacaaggatcgagaccatgacaggggccgggagaagagatggcaagaagagagccgaccagggtgtggccgacaatgactgggagagagagggacttcagagat
gacaggatcaaggaggaggagaagaaggaaaggaggaagtag