



Reference: PTXBC009684

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MCT5; MCT6; monocarboxylate transporter 6; monocarboxylate transporter 5; solute carrier family 16 (monocarboxylate transporter), member 5; solute carrier family 16 (monocarboxylic acid transporters), member 5; solute carrier family 16, member 5 (monocarb		

For research use only.

GENE DESCRIPTION	solute carrier family 16, member 5 (monocarboxylic acid transporter 6)
GENE ID	9121
ORIGIN SPECIES	Human
ACCESSIONS	BC009684
VECTOR	pDONR223
NCBI SYMBOL	SLC16A5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]