

DNAJC19-DnaJ (Hsp40) homolog, subfamily C, member 19 Gene

Reference: PTXBC009702

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAM18; TIM14; TIMM14; mitochondrial import inner membrane translocase subunit TIM14; DnaJ (Hsp40) homolog, subfamily C, member 19; DnaJ-like protein subfamily C member 19; homolog of yeast TIM14; DnaJ heat shock protein family (Hsp40) member C19		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 19
GENE ID	131118
ORIGIN SPECIES	Human
ACCESSIONS	BC009702
VECTOR	pDONR223
NCBI SYMBOL	DNAJC19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagtacagtggtagcagttggactgaccattgctgctgcaggatttgcaggcgttacgttttcaagccatgaagcatatggagcctcaagtaaaacaagttttcaagcctacaaaatctgccttcagtggtggctattatagaggtgggttgaacccaaaatgacaaaacggg
aagcagcattatatactaggtgtaagccctactgccaataaagggaataaagagatgctcatcgacgaattatgcttttaaatcatcctgacaaggaggtctccttatatagcagccaaaatcaatgaagctaaagatttactagaaggtcaagctaaaaatga