

SEC11C-SEC11 homolog C (*S. cerevisiae*) Gene

Reference: PTXBC009703

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	signal peptidase complex catalytic subunit SEC11C; SEC11L3; SPC21; SPCS4C; SEC11 homolog C; SEC11-like 3; SEC11-like protein 3; SPase 21 kDa subunit; microsomal signal peptidase 21 kDa subunit; signal peptidase complex 21; SEC11 homolog C, signal peptidas		

For research use only.

GENE DESCRIPTION	SEC11 homolog C (<i>S. cerevisiae</i>)
GENE ID	90701
ORIGIN SPECIES	Human
ACCESSIONS	BC009703
VECTOR	pDONR223
NCBI SYMBOL	SEC11C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcgtgcggcgccgtggggctcatctccccgcgtccgctggatctctggggacctgaagaagatgaacaagcgccagctctattaccaggttttaacttcgccatgatcgtgtctctgcactcatgatatgaaaggcttgatcgtgtcacaggcagtgagagcccatc
gtgggtgctgagtgagcagtgatggagccggccttcacagaggagacctctctcacaatttcgggaagaccaatcagagctggtgaaatagttgttttaagggtgaaggacgagacattccaatagttcacagagtaataaagttcatgaaaagataatggagacataaatt
tctgactaaaggagataataatgaagttgatagaggtgttacaaagaaggcagaactggctgaaaagaaggacgtgggtgggaagagcaagaggggttttaccatatgttggtatggtcaccataataatgaatgactatccaaatcaagtatgctctttggctgtaatgggtgcat
atgtgttactaaaacgtgaatcctaa