



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TSPAN5-tetraspanin 5 Gene

Reference: PTXBC009704

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET-4; NET4; TM4SF9; TSPAN-5; tetraspanin-5; tetraspan NET-4; tetraspan TM4SF; transmembrane 4 superfamily member 9; transmembrane 4 superfamily, member 8; tetraspanin 5		

For research use only.

GENE DESCRIPTION	tetraspanin 5
GENE ID	10098
ORIGIN SPECIES	Human
ACCESSIONS	BC009704
VECTOR	pDONR223
NCBI SYMBOL	TSPAN5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcgcggagacacacatacgaaggtcctgaagtcagttgttcacataaactcatttttggcttcaagtctatatttttgggtttgggaataacattcttggaattggagctgtggcgatggaatgaaaaaggagttctccaacatctctccatcacacgtctggcggttgaccagtttggctcttcctgtgtgtggggaggtgatgttcattttgggatttcagggtgcatggagcgctacgggaaaaacattctctcacaagttttttctgtgtctctgggaattatttctctctggagctcactcggcgagttctacattgttttcaaagactggaatcaaaagaccagctgtatttcttataaacaacacatcacagatcctgggatgacattgatttgcacaaacctatagacttcaccaggaaatttgcagatgctgtgtgggctttggagagctgatgattggaacctaataatttactcaattgcacagatttccaatgcaagtcgagagcgaattggcggttcattctctctgtgcactaaagatcccgacagaagattgtcatcacactcagttgtgctatgatgcagtcgacaaacagaagtgaccacagattgtaacttacacgaagagctgtgtgcccgattgagaagttgtgtcaggacaafttaaccatggtctgtgtatttcatagcgattgctgctgcagatattgggatagcttggcccaaaatgttcgcacatacgaagctctcagtcgcagctgctgtag

More info: Please contact us for any question: contact@proteogenix.fr